^{en} 23G Cannula (1) ^{da} Kanyla 23G (1) ^{nl} Canule, 23 G (1)
^{fr} Canule 23G (1) ^{de} 23G Kanüle (1) ^{it} Cannula da 23 G (1)
^{no} 23G kanyla (1) ^{pt} Cãnula de calibre 23 G (1) ^{es} Cãnula de 23G (1) ^{sv} 23G Kanyl (1) ^{ar} قنية 23
^{bg} Канюла 23G (1) ^{cs} 23G kanyla (1) ^{fi} 23 G:n kanyyli (1)
^{el} Σωληνίσκος 23G (1) ^{he} צינורית 23G (1) ^{hu} 23G kanül (1)
^{pl} Kaniula 23G (1) ^{ro} Canulã de 23G (1) ^{ru} Канюля 23G (1) ^{sk} Kanyla 23G (1) ^{tr} 23G Kanül (1) th 23G ฟิล์นาล์ (1)
^{uk} Канюля 23G (1) ^{lv} 23G kanüül (1) ^{lt} 23G kaniulė (1)
^{lv} 23G kanula (1)

^{en} MANI Ophthalmic Knife (1mm Slit Angle Knife) (1)
^{da} MANI oftalmisk kniv (kniv med 1 mm spaltevinkel) (1)
^{nl} MANI Ophthalmic Knife (Gehoek slit-mesje, 1 mm) (1)
^{fr} Scalpel ophthalmique MANI (scalpel biseauté de 1 mm) (1)
^{de} MANI Ophthalmisches Messer (1 mm Schlitz-Winkel-Messer) (1)
^{pt} Bisturi a punta angolata da 1 mm MANI Ophthalmic Knife (1)
^{es} MANI oftalmisk kniv (1 mm spaltevinkelkniv) (1)
^{sv} Bisturi oftálmico MANI (bisturi de corte em ângulo de 1 mm) (1)
^{cs} Cuchillete oftálmico MANI (cuchillete angulado de hendidura de 1 mm) (1)
^{tr} MANI oftalmisk kniv (1 mm slitsvinkelkniv) (1) (1) (م 1 Slit Angle مشرط عين)
^{bg} Офталмологичен скалпел MANI (ъглов слит-скалпел 1 mm) (1)
^{en} Oftalmologický nůž MANI (1mm štrbinový zahnutý nůž) (1)
^{fi} MANI-silmäveitsi (1 mm:n kulmaveitsi) (1)
^{el} Οφθαλμική Λεπίδα MANI (Γωνιακή Λεπίδα Λεπτής σχισμής 1 mm) (1)
^{he} סכין מני (סכין לחתוך בזווית בעובי של 1 מ"מ) (1)
^{hu} MANI szemészeti szike (1 mm-es résszögű szike) (1)
^{pl} Skalpel okulistyczny MANI (skalpel kątowy 1 mm) (1)
^{ro} Bisturiu oftalmic MANI (bisturiu unghiular cu pas dublu de 1 mm) (1)
^{ru} Офтальмологический нож MANI (расплаиватель с двойной заточкой 1 мм) (1)
^{sk} Oftalmický nůž MANI (1 mm štrbinový zahnutý nůž) (1)
^{tr} MANI Oftalmik Bıçak (1 mm Slit Açılı Bıçak) (1)
th MANI ฟิล์นาล์ดัดปลาย ฟิล์นาล์ (1) (ฟิล์นาล์ดัดปลาย ฟิล์นาล์) (1)
^{uk} Офтальмологічний ніж MANI (зігнутий для тунельного розрізу 1 мм) (1)
^{en} MANI silmanuga (1 mm, sisseõikeks, kaldteraga) (1)
^{pl} MANI oftalmologinis peilis (1 mm plyšio kampo peilis) (1)
^{lv} MANI oftalmoloģijas nazis (1 mm spraugas leņķveida nazis) (1)

PRESERFLO™ MicroShunt GLAUCOMA DRAINAGE SYSTEM (en)
Instructions for Use

INTRODUCTION

The PRESERFLO™ MicroShunt product consists of the PRESERFLO MicroShunt and a 3-20mm stainless steel Scleral Marker. The MicroShunt is an implantable glaucoma drainage device made of an extremely flexible SIBS [poly (Styrene-block-IsoButylene-block-Styrene)] polymer with a tube of 350 µm outer diameter and a lumen of 70 µm. It has triangular fins that prevent migration of the tube into the AC. The device is designed to be implanted under the subconjunctival/Tenon space.

1. Make sure to check the expiration date on each item.

INTENDED USE/ PURPOSE

The PRESERFLO MicroShunt employs a tube to create a conduit for the flow of aqueous humor from the anterior chamber of the eye to a bleb formed under the conjunctiva and Tenon's capsule; the front/ proximal end of the tube extends into the anterior chamber while the back/ distal end terminates in the bleb. The PRESERFLO MicroShunt reduces IOP by physically shunting aqueous from the high pressure anterior chamber to the lower pressure bleb.

INTENDED USERS

The intended users of the PRESERFLO MicroShunt are ophthalmologists/ophthalmic surgeons specializing in the treatment of glaucoma (including surgeons specializing in anterior segment and cataract surgery) who have been trained to use the device.

INDICATIONS FOR USE

The PRESERFLO™ MicroShunt Glaucoma Drainage System is intended for reduction of intraocular pressure in eyes of patients with primary open angle glaucoma where IOP remains uncontrollable while on maximum tolerated medical therapy and/or where glaucoma progression warrants surgery.

INTENDED PATIENT/ TARGET GROUPS

The target patient population is adult patients with primary open-angle glaucoma where IOP remains uncontrolled while on maximum tolerated medical therapy and/or where glaucoma progression warrants surgery.

CONTRAINDICATIONS

The implantation of the PRESERFLO™ MicroShunt is contraindicated under the following circumstances and conditions: ANGLE CLOSURE GLAUCOMA; PRESENCE OF CONJUNCTIVAL SCARRING, PREVIOUS INCISIONAL OPHTHALMIC SURGERY INVOLVING THE CONJUNCTIVA OR OTHER CONJUNCTIVAL PATHOLOGIES (E.G., THIN CONJUNCTIVA, PTERYGIUM) IN THE TARGET QUADRANT; ACTIVE IRIS NEOVASCULARIZATION; ACTIVE INFLAMMATION (E.G., BLEPHARITIS, CONJUNCTIVITIS, SCLERITIS, KERATITIS, UVEITIS); VITREOUS IN THE ANTERIOR CHAMBER; PRESENCE OF AN ANTERIOR CHAMBER INTRAOCULAR LENS (ACIOL); INTRAOCULAR SILICONE OIL.

WARNINGS

- Rx only: This device is restricted to sale by, or on the order of, a physician.
- For one-time use only. Do not reuse or re-sterilize. Reuse, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in serious patient injury, illness, blindness or death. Reuse, or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, blindness, or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.
- Long term effects of Mitomycin C (MMC) with the use of this device have not been evaluated. Necessary precautions and interventions on the use of MMC are highly recommended.
- The PRESERFLO™ MicroShunt should not be subjected to direct contact with petrolatum-based (i.e., petroleum jelly) materials (e.g., ointments, dispersions, etc.).
- The effects of cutting or modifying the PRESERFLO™ MicroShunt have not been evaluated.
- Viscoelastics have not been tested with the PRESERFLO™ MicroShunt. However, in an emergency when all other therapies have failed, the use of hydroxypropyl methylcellulose (HPMC) may be an option. Use of HPMC should be a last resort to correct a flat chamber with the PRESERFLO™ MicroShunt and may risk loss of flow through the device for one or more weeks after use necessitating close or more frequent observation of IOP.

PRECAUTIONS

- The safety and effectiveness of the PRESERFLO™ MicroShunt has not been established in patients with the following circumstances or conditions:
 - o Eyes with chronic inflammation.
 - o Eyes with congenital and infantile glaucoma.
 - o Eyes with neovascular glaucoma.
 - o Eyes with uveitic glaucoma.
 - o Eyes with pseudoexfoliative or pigmentary glaucoma.
 - o Eyes with other secondary open angle glaucomas.
 - o Eyes that have undergone previous incisional glaucoma surgery or cilioabative procedures.
 - o Implantation with concomitant cataract surgery with IOL implantation.
- Avoid use of toothed forceps to handle device. McPherson type forceps are recommended.
- The patient's IOP should be monitored postoperatively. If IOP is not adequately maintained after surgery, appropriate additional therapy to maintain IOP should be considered.
- The safety and effectiveness of the use of more than a single PRESERFLO™ MicroShunt has not been established.
- If the PRESERFLO™ MicroShunt appears deformed, folded and/or distorted, do not use.
- The creation of 2 or more scleral tunnels in an attempt to implant the MicroShunt may cause leakage of aqueous humor and increase the risk of hypotony, if the additional scleral tunnels did not self-seal or sutured.

- When managing patients with a shallow or flat anterior chamber, while not a likely occurrence, there is a possibility that pharmacologic dilation of the pupil in these patients may cause the iris to obstruct the MicroShunt. For these patients, more frequent monitoring of IOP following dilation is recommended.

POTENTIAL COMPLICATIONS/ADVERSE EVENTS

The complications during and after surgery may include:

Glaucoma progression not controlled, increase in cup-to-disc ratio (C/D), anesthesia related complications, difficulty in inserting the PFMS, failure to implant the device, device malfunction, device repositioning, extended surgical procedure, tube migration out of anterior chamber, flat anterior chamber, shallow anterior chamber, excessive bleeding in anterior chamber or incision site or eye, PFMS touches cornea or iris, intraocular pressure too high or low, strabismus, viscoelastic used in anterior chamber, choroidal effusion or hemorrhage, retinal complications (retinal detachment, proliferative retinopathy), hyphema, hypotony or hypotony maculopathy, phthisis bulbi, endophthalmitis, tube erosion through conjunctiva, tube obstruction partial or complete (block by iris or vitreous or fibrin or debris), iridodialysis, iris incarceration, uveitis, iritis, diplopia, aqueous misdirection, corneal complications (abrasion, edema, ulceration, infection, decompensation, bullous keratopathy, endothelial cell loss, Descemet striae, keratitis), partial or complete vision loss, blurry vision, bleb leak, blebitis, cystic bleb, bleb failure, pupillary block, ptosis, macular edema, prolonged inflammation, use of glaucoma medications, ocular pain, conjunctival complications (buttonhole, dehiscence, dissection, hemorrhage, hyperemia, scar, tear, ulcer), iris adhesions/synechiae, iris abrasion, cataract development or progression, posterior capsule opacity, explantation of the PFMS, encapsulation reaction, foreign body sensation, fibrin in anterior chamber, visual field damage, unplanned glaucoma-related surgical re-intervention, optic disc hemorrhage, globe perforation, headache, vitreous hemorrhage, wound leak, and suture related complications.

STERILIZATION

All PRESERFLO™ MicroShunt devices and accessories are sterile upon receipt. The PRESERFLO™ MicroShunt, scleral marker, and Slit Angle Knife are sterilized by ETO sterilization per ISO 11135 guidelines and the marker pen, needle and cannula, are sterilized by gamma radiation. Contents will remain sterile unless package is opened or damaged. Should the device or any accessory be inadvertently rendered non-sterile, or should the package be opened and/or damaged, discard the device. **DO NOT RESTERILIZE.**

STORAGE

Store at controlled room temperature (20±5°C) and avoid excessive heat, excessive humidity, direct sunlight, and water.

INSTRUCTIONS FOR USE

User must be properly trained on the use of the device prior to use.

Surgery:

- Anesthesia:** The type of anesthesia to be administered during the procedure is at the surgeon's discretion.
 - Preparation of PRESERFLO™ MicroShunt:** a) Carefully examine the pouch containing the PRESERFLO™ MicroShunt. If pouch is damaged, do not use the device. b) Remove PRESERFLO™ MicroShunt from sterile packaging onto sterile field. c) Wet the PRESERFLO™ MicroShunt using a solution of Balanced Salt Solution (BSS).
 - Planning of Implant Position:** Insertion of the PRESERFLO™ MicroShunt is at the surgeon's discretion, either in the superonasal or superotemporal quadrant of the eye, avoiding the rectus muscles.
 - Corneal Traction and Conjunctival Flap:** Perform corneal traction with suture at surgeon's discretion. In the planned surgical area, make a 6-8mm peritomy incision, then dissect a subconjunctival/subTenon pouch (fornix based) over a circumference of 90°-120° from near the limbus to at least 8mm deep. Assure there is no active bleeding in the area of the planned PRESERFLO™ MicroShunt entry. Avoid excessive use of cautery, which could increase the risk of scarring.
 - Application of MMC:** At the surgeon's discretion, apply sponges saturated with MMC on the sclera under the subconjunctival/subTenon pouch, making sure they do not come into contact with the conjunctival wound edge. Remove the sponges and rinse with generous amounts of BSS (>20mL).
 - Preparation of Anterior Chamber Entrance Site:** Ink the 3-20mm Marker (Scleral Marker) points with the Marker Pen - Gentian Violet and mark a point 3mm from the limbus.
 - If using the **MANI Double Step-Angled 1.0mm Knife:** At the marked point, create a scleral tunnel into the anterior chamber (AC) entering just above the iris plane (level of trabecular meshwork). As a guide externally, enter the AC just below the middle border of the surgical limbus (posterior border of blue/gray zone, which should overlie internally with Schwalbe's line), if present. The scleral tunnel is complete when a scleral pocket that is 2mm long is created at the distal end of the scleral tunnel, where the fin portion of the PRESERFLO™ MicroShunt is tucked.
 - If using the **MANI Slit Angled 1.0mm Knife:** At the marked point, create a shallow scleral pocket that is 2mm long, where the fin portion of the PRESERFLO™ MicroShunt is tucked. Insert a 25G needle into the scleral pocket and create a scleral tunnel into the AC entering just above the iris plane (level of trabecular meshwork). As a guide externally, enter the AC just below the middle border of the surgical limbus (posterior border of blue/gray zone, which should overlie internally with Schwalbe's line), if present.
 - Insertion of PRESERFLO™ MicroShunt into Anterior Chamber: Using a non-toothed forceps,** position the PRESERFLO™ MicroShunt tube with the bevel facing towards the cornea close to the scleral pocket/scleral tunnel while keeping the fin face parallel to the scleral surface. Gently insert the PRESERFLO™ MicroShunt into the scleral pocket/scleral tunnel until the fin is tucked under the scleral pocket. Use of excessive force will result in temporary bending of the PFMS and increase the difficulty of entrance. If not accessible, another scleral pocket/scleral tunnel can be made approximately 1mm to either side of the original tract. If the PRESERFLO™ MicroShunt is damaged during the procedure, use a new PRESERFLO™ MicroShunt and use a new scleral pocket/scleral tunnel. After successful insertion of the PRESERFLO™ MicroShunt, ensure that the proximal end of the tube in the AC is not in contact with the cornea or iris, and that there is flow of aqueous from the distal end. Then tuck the end of the PRESERFLO™ MicroShunt tube under the subconjunctival/subTenon's flap on the scleral surface, making sure it is straight and not caught in Tenon's capsule. If flow is not visible, consider one or more of the following actions: applying gentle pressure on the eye (sclera, limbus, cornea), filling the anterior chamber with BSS through a paracentesis, or using a 23G Thin wall 8mm bend cannula and inject BSS into the distal end of PRESERFLO™ MicroShunt to prime the tube. If there is no flow established, reposition or replace the device.
 - Verify position of the MicroShunt in the anterior chamber angle:** It is recommended that intraoperative gonioscopy be performed to assess the position of the proximal tip of the PRESERFLO™ MicroShunt in the anterior chamber to ensure that there is sufficient clearance from the cornea or iris.
 - Conjunctival Closure:** Reposition the Tenon's and conjunctiva to the limbus and suture using sutures with a long history of successful use by glaucoma specialists (i.e., nylon, vicryl). A moistened fluorescein strip is used to check for a conjunctival leak, which may be closed with additional sutures. Verify the following: presence of the proximal end of the PRESERFLO™ MicroShunt in the anterior chamber; distal end of the PRESERFLO™ MicroShunt is straight and not caught in Tenon's capsule. Pull corneal traction suture out if it was used. Apply antibiotic and steroidal medication postoperatively as appropriate.
 - After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy. Follow the MMC instructions for use for proper disposal of the MMC.
- Post operatively:**
- Monitor the intraocular pressure at each follow-up to determine if the PRESERFLO™ MicroShunt is patent.
 - If the PRESERFLO™ MicroShunt is repositioned, removed and/or replaced with another type of device based on the device performance, the conjunctiva should be cut at the limbus in a similar manner to the original procedure. The device should be exposed, and the repositioning or removal conducted. If the device is removed, verify that there is no aqueous leakage from the track. If there is, suture the track to seal it. Conjunctival closure should follow the original procedure.

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PRESERFLO MicroShunt (MicroShunt & 3 mm Scleral Marker)

Manufacturer of the PRESERFLO™ MicroShunt

INNFOCUS

A Santen Company
12415 S.W. 136 Avenue
Miami, FL, 33186
USA
Phone: +1(305) 378-2651
Fax: +1(305) 378-2652
www.innfocusinc.com

EC REP

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
www.emergobyul.com

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UK Responsible Person

Santen UK Ltd
Salisbury Hall, St Albans, AL2 1BU,
United Kingdom

NOTE: InnFocus is the Manufacturer of the PRESERFLO™ MicroShunt and the assembler of the Procedure Pack.

Date of Issue: November 21, 2024

CAW-4147-05

For any medical inquires or safety reports contact Santen Medical Information Call Center at +1 (855) 772-6836 or by email at medicalinformation@santen.com or contact the local Santen Office at www.santen.eu. Any serious incident that has occurred in relation to the device should be reported to Santen and the competent authority of the Member State in which the user and/or patient is established.

Order	Language	ISO Abbreviation
1	English	en
2	Danish	da
3	Dutch	nl
4	French	fr
5	German	de
6	Italian	it
7	Norwegian	no
8	Portuguese	pt
9	Spanish	es
10	Swedish	sv
11	Arabic	ar
12	Bulgarian	bg
13	Czech	cs
14	Finnish	fi
15	Greek	el
16	Hebrew	he
17	Hungarian	hu
18	Polish	pl
19	Romanian	ro
20	Russian	ru
21	Slovak	sk
22	Turkish	tr
23	Armenian	hy
24	Ukrainian	uk
25	Estonian	et
26	Lithuanian	lt
27	Latvian	lv